

gandio™



Treatment - Simple, Fast and Safe

A daily use of 5 minutes for users with short hair and 8 minutes for those with longer hair will bring a significant improvement within several months. It is recommended to use 1-2 times daily, however, for faster and better results use it twice a day for the first two months.

Once major regrowth potential has been achieved you may reduce the treatment to 3-4 times per week.

1. Lift the Applicator from the Base.
2. Lightly press the Applicator's button. After a quick self-test, the Applicator will emit a startup tune and start vibrating.
3. For effective treatment, roll the Applicator back and forth over the scalp area requiring treatment.
4. Profile II - 5 minutes treatment - recommended for short hair
5. Profile III (8 Minutes) is intended for people with long or dense hair. For such people, it is recommended to part the hair and roll back and forth across the parting line for half a minute, then move on to another parting line.
6. If you find the intensity and sensitivity too strong, you can use profile I to let the scalp get use to the treatment.
7. The Base and Applicator should be cleaned with a damp cloth and then dried. The disks and the area between them can be cleaned with a brush. For sanitary purposes, roll the disks on an alcohol pad, and dry before use.